



Glow 70mg PASS

Tested for: Prism Labs
prismlabs.store

COA #
PPA-COA-2026-575805

Lot Number
575805

Accession #
PPA-2026-575805

Concentration
70mg

Sample Matrix
Lyophilized

Method
Full QC Panel

Analysis Date
08/03/2026

Appearance
Good

Volume
3mL

Received
07/27/2026



Scan to verify authenticity

<https://www.proofpointanalytical.com/coa-lookup/?code=YXF9OUVH>

Identity

Glow-70

Purity

99.25%

Blend Purity & Quant (HPLC)

Analyte	Specification	Result	Unit	Status
Peptide Purity (HPLC)	>= 95.0%	99.25%	%	PASS
Net Blend Peptide Content	Report Only	72.53	mg	N/A
- GHK-Cu		52.05	mg	
- BPC-157		10.48	mg	
- TB-4		10.00	mg	
Identity (ID)	GHK-Cu + BPC-157 + TB-4	Confirmed	-	PASS

Authorized Review
Lab Director

Issued: 08/03/2026 | Access Code: YXF9OUVH



Glow 70mg PASS

COA #	PPA-COA-2026-575805	Lot Number	575805
Accession #	PPA-2026-575805	Analysis Date	08/03/2026
Access Code	YXF9OUVH	Verify	https://www.proofpointanalytical.com/coa-lookup/?code=YXF9OUVH

HPLC Conformity Testing Results

Sample	Purity	NPC	ID Result	Status
Dedicated V0	99.13%	72.55 mg	Confirmed	PASS
Conformity V1	99.25%	72.53 mg	Confirmed	PASS
Mean	99.19%	72.54 mg	-	-
Std Dev	0.0600%	0.0100 mg	-	-

Heavy Metals Analysis (ICP-MS)

Test	Specification	Result	Status
Arsenic (As)	NMT 1.5 ppm	Not Detected	PASS
Cadmium (Cd)	NMT 0.5 ppm	Not Detected	PASS
Chromium (Cr)	NMT 10 ppm	Not Detected	PASS
Mercury (Hg)	NMT 1.5 ppm	Not Detected	PASS
Lead (Pb)	NMT 1 ppm	Not Detected	PASS

Sterility Testing (PCR)

Test	Specification	Result	Status
Sterility (PCR)	No Growth	No Growth	PASS



Glow 70mg PASS

Endotoxin Testing (USP <85>)

Test	Specification	Result	Status
Endotoxin (USP <85>)	Report Result	NMT 0.05 EU/mL	Reported

About this result: Endotoxin is reported as a quantitative screening value. Acceptable limits vary by product type and matrix, so no universal pass/fail threshold applies across product types.

Notes & Methodology

- 1. Methods: Purity & Quant (HPLC), identity confirmation, heavy metals screening, sterility PCR, and endotoxin review as applicable.
- 2. Peptide purity determined by RP-HPLC area normalization at 214 nm. Value represents the percentage of the target peptide relative to all peptide-related peaks.
- 3. Identity confirmed by chromatographic retention time comparison and analytical review against expected product identity.
- 4. Elemental impurities analyzed by ICP-MS screening methodology using internal quality limits for research-use materials.
- 5. Results relate only to the item tested and are intended for research-use-only materials.

Access Code: YXF9OUVH

Issued: 08/03/2026



Verify COA